ASPINA - DO M

NÁVOD NA INŠTALÁCIU, OBSLUHU A ÚDRŽBU
INSTALLATION, OPERATION AND MAINTENANCE MANUAL
ИНСТРУКЦИЯ ПО УСТАНОВКЕ, ОБСЛУЖИВАНИЮ И УХОДУ
DIE ANLEITUNG ZUR INSTALLATION, BEDIENUNG UND WARTUNG
NOTICE D'INSTALLATION, DE MISE EN MARCHÉ ET D'ENTRETIEN
INSTRUKCJA INSTALACJI, OBSŁUGI ORAZ NAPRAWY BIEŻĄCEJ

Ⓚ

ⒼⒷ

ⓇⓊ

Ⓓ

Ⓕ

Ⓟ

CE 0434

ekom[®]

**MEDICINÁLNY KOMPRESOR
MEDICAL COMPRESSOR
МЕДИЦИНСКИЙ КОМПРЕССОР
MEDIZINALKOMPRESSOR
LE COMPRESSEUR À USAGE MÉDICAL**

**VÝROBCA
MANUFACTURER
ПРОИЗВОДИТЕЛЬ
HERSTELLER FABRICANT
PRODUCENT**

**DÁTUM POSLEDNEJ REVÍZIE
DATE OF LAST REVISION
ДАТА ПОСЛЕДНЕГО ПЕРЕСМОТРА
DATUM DER LETZTEN ÜBERARBEITUNG
DATE DE MISE A JOUR
DATA WYDANIA**

Aspina -DO M

EKOM spol. s r. o.
Priemyselná 5031/18
SK-921 01 Piešťany
Slovak Republic
tel.: +421 33 7967255
fax: +421 33 7967223

www.ekom.sk
email: ekom@ekom.sk

01/2012

SK	<u>OBSAH..... 2</u>
GB	<u>CONTENTS..... 13</u>
RU	<u>СОДЕРЖАНИЕ..... 24</u>
D	<u>INHALT..... 37</u>
F	<u>SOMMAIRE..... 48</u>
PL	<u>SPIS TREŚCI..... 60</u>

OBSAH

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 1.2. General warnings 14

 1.3. General safety warnings 14

 1.4. Safety warnings regarding the protection against electric current 14

 1.5. Alert notices and symbols 14

2. PRODUCT INFORMATION 15

 2.1. Intended Use 15

 2.2. Product Description 16

3. TECHNICAL DATA 16

4. FUNCTION DESCRIPTION 17

 4.1. Function Description 17

 4.2. Detailed Description of the Suction Part Function 17

ASSEMBLY 17

5. STORAGE AND TRANSPORT CONDITIONS 17

6. INSTALLATION AND INITIAL START-UP 17

 6.1. Environmental Conditions 18

 6.2. Electric Connection 18

 6.3. Initial Start-up 18

USE 18

7. OPERATION 18

 7.1. Switching on the Dental Aspirator 18

 7.2. Use of antifoaming tablets 19

8. MAINTENANCE INTERVALS - USER / TECHNICIAN 19

9. MAINTENANCE, CLEANING AND DISINFECTION 20

 9.1. Inlet Sieve Cleaning 20

 9.2. Disinfection of Tubings and Separation Vessel 20

 9.3. Cleaning and disinfection of the exterior surfaces of the product 20

 9.4. Replacement of the Output Pre-filter 20

 9.5. Replacement of Output Filter 21

 9.6. Sterilisation Process for Cannulas (pursuant to STN EN ISO 17664) 21

10. PUTTING OUT OF OPERATION 22

DISPOSAL 22

11. DISPOSAL OF APPLIANCE 22

TROUBLESHOOTING 22

12. INFORMATION ON REPAIR SERVICE 22

13. SOLVING COMMON PROBLEMS 22

14. WIRING DIAGRAMS 71

15. FUNCTION DIAGRAM 73

ENCLOSURE NO. 1 74

SCOPE OF DELIVERY 75

GUARANTEE 76

IMPORTANT INFORMATION**1. WARNINGS****1.1. CE Marking**

Products marked with **CE** mark of compliance meet safety guidelines of European Union (93/42/EEC).

1.2. General warnings

- The installation, operation and maintenance manual is an integral part of the appliance. It is necessary to always keep this document close to the appliance. Strict observance of this manual is a prerequisite for the correct operation of the appliance.
- The safety of operating personnel and failure-free operation of the appliance are ensured only when using the original components of the appliance. It is possible to use accessories and spare parts mentioned in the technical documentation or clearly permitted by the manufacturer.
- When used with non authorized accessories or consumable material, the manufacturer cannot assume responsibility for the safe operation and functionality of the device.
- The Guarantee does not cover damages that originate due to the use of non authorized accessories or consumable material other than those recommended by the manufacturer.
- The manufacturer assumes responsibility regarding safety, reliability and function only if
 - The installation, calibration, amendments, extensions and repairs are made by the manufacturer or his representative or a service organization authorized by the manufacturer,
 - The appliance is used in accordance with the installation, operation and maintenance manual.
- The installation, operation and maintenance manual, at the time of printing, corresponds to the design of the appliance and its state according to the relevant safety and technical standards. The manufacturer reserves all copyrights for the given wiring schemes, methods and names.

1.3. General safety warnings

The manufacturer developed and constructed the appliance so that damage would not occur when the appliance is used for its intended purpose. The manufacturer considers it his obligation to describe the following safety measures in order to avoid further damages.



- When operating the appliance, it is necessary to observe laws and regional regulations valid in the place of usage. In order to ensure safe course of works, the operator and user are responsible for the observation of regulations.
- The original packaging should be kept for the possible return of the unit. Only original packaging guarantees an optimal protection of the appliance during transportation. If it would be necessary to return the appliance during warranty period, the manufacturer is not responsible for damages caused by incorrect packing.
- It is necessary that the user ensures the appliance is safe to use prior to usage.
- The user must familiarize himself with the correct operation of appliance.
- If an undesirable event occurs in the operation of appliance, the user is obliged to immediately inform his supplier to this event.
- This product is not intended for use in areas with the risk of explosion.

1.4. Safety warnings regarding the protection against electric current

- The appliance must only be connected to an appropriate power source that has correct grounding.
- Prior the connecting the compressor, verify whether the mains voltage and frequency specified on the apparatus are in accordance with the local supply.
- Prior to putting into operation, check for possible damages on the appliance and the air connectors. Damaged cables and sockets/plugs must be replaced immediately.
- In the case of a dangerous situation or a technical failure, immediately disconnect the appliance from mains supply.
- During all repairs and maintenance:
 - ensure that the mains plug is removed from the power socket
 - pressure pipes must be air vented
 - pressure must be released from pressure tank.
- This appliance can only be installed only by a qualified expert.

1.5. Alert notices and symbols

For your information, the symbols below are used in the installation, operation and maintenance manual, on packaging materials and on the product :

	Information or instructions to prevent any injury to health or material damage to the compressor.
	Alert against dangerous electric voltage.
	Consult instructions for use
	CE – marking
	Attention! Hot surface.
	Handling mark on packaging – Fragile, handle with care
	Handling mark on packaging – This way up (Vertical position of freight)
	Handling mark on packaging – Protect against humidity
	Handling mark on packaging – Temperature of storage and transport
	Handling mark on packaging – Limited stacking
	Mark on packaging – Recyclable material
	Connection of protective ground wire
	Connector for equipotential connection
	Fuse
	Alternating current
	Danger of biological hazard

2. PRODUCT INFORMATION

2.1. Intended Use

The dental aspirator ASPINA DO M is a mobile equipment extending the possibilities of procedures provided for by a physician and ensuring the increase of ergonomcy at his work. It is suitable for the dental units not equipped with the suction apparatus and a separator. It is intended for the dental practice needs for aspiration, separation and entrapping of the waste into the built-in separation vessel. With regard to its mobility it enables a simple transport between the workplaces and can be used above all where the connection to a sewerage system is impossible.

- The dental aspirator is designed for the operation in dry, ventilated rooms, with the ambient temperature ranging between +5°C up to +40°C, and relative humidity not exceeding the value of 70%.

- The dental aspirator may not be exposed to rain. The apparatus may not be operated in the moist or wet environs. In addition, avoid the use near gases or combustible liquids.
- Other use, or the use exceeding this scope, can not be regarded as the intended use. The manufacturer is not responsible for the resulting damages. The risk will be born exclusively by the operator / user.

2.2. Product Description

The mobile dental aspirator ASPINA DO M is built-up on a movable bogie on which is situated a box, lined by a noise damping material. Inside the box there is situated a fan cooled suction aggregate (9) with the electric distribution and the waste entrapping separation vessel (11). In the lower part – under the box – there is situated a silencer with an output filter (14) and a pre-filter (15) ensuring the air filtration from the suction aggregate. In the upper – narrowed part of the aspirator - there is situated a holder for suction tubings (2) equipped with mouthpieces (1), separation automatics and terminal box with fuses. On the lateral part there is situated a main switch (5), over which are placed the indicators for the network (3) and for the state of filling of the separation vessel (4).

3. TECHNICAL DATA

	DO M
Rated voltage / frequency V / Hz	230 / 50; 230 / 60*; 110 / 60*
Rated input of the aspirator VA	520
maximum flowrate : l/min	1100
maximum underpressure kPa	12
Sound level dB(A)	≤48
Mode of operation of the aspirator	permanent S 1
Dimensions of the aspirator mm	565x350x860
Dimensions of the aspirator with carton (h x w x d) mm	950x560x595
Weight of the aspirator kg	36
Weight of the aspirator with carton kg	40
Model according to EN 60 601-1 (IEC 601-1) STN	apparatus type B, class I.

(*) The rated voltage and frequency can be change by special order

Climatic conditions for the storage and transport

Temperature –25°C up to +55°C, 24 h up to +70°C

Relative air humidity 10% up to 90 % (without condensation)

Climatic conditions for the operation

Temperature +5°C up to 40°C

Relative air humidity up to +70%

4. FUNCTION DESCRIPTION

4.1. Function Description

After switching-on the main switch (5) into the position „I“ the network indicator (3) goes on. Taking the suction tubing (6) off the holder (2) will actuate the suction aggregate (9) and the underpressure occurs at the suction mouthpiece (1). After repeated putting the suction tubing into the holder, the suction aggregate turns off. When the separation vessel (11) is filled with waste products, the suction aggregate turns off and the indicator for the separation vessel (4) filling goes on. Then it is necessary to return the suction tubing back into the holder and to empty the separation vessel. During a longer work, mainly with the desalivating mouthpiece, the box temperature may increase; the cooling fan is then automatically actuated. The fan turns off automatically, when the box temperature drops

4.2. Detailed Description of the Suction Part Function

The underpressure air flows together with sucked waste products from the oral cavity through the tubing system from the suction mouthpiece (1) at first through the inlet sieve (7), wherein the solid impurities are entrapped. The underpressure air, together with the sucked waste products, free of solid impurities greater than 2 mm, then flows into the separation vessel (11) wherein it is separated from the underpressure air and entrapped into the separation vessel.

The underpressure air flows afterwards into the suction aggregate (9), from where it is blown off through the silencer. Therein the output air passes through the output pre-filter (15) and the bacteriologic output filter (14). After passing the filters, the air free of impurities is blown off into a free space under the dental aspirator.

ASSEMBLY

5. STORAGE AND TRANSPORT CONDITIONS

The dental aspirator is delivered ex-work packed in a transport carton, preventing the apparatus from being damaged during the transport.

**For the transport, always use an original product packing, if possible.
Transport the dental aspirator in an upright position.**

Protect the aspirator during transport against moisture, impurity or extreme temperatures. The aspirator in an original packing can be stored in warm, dry and dustfree rooms.

Keep the packing material for possible future return shipments. If this is not possible, dispose the packing material in an environmental-friendly way. The transport carton can be disposed off as the old paper.

**Mobile dental aspirator may only be transported with emptied separation vessel.
Always empty the separation vessel contents before any transport.**

6. INSTALLATION AND INITIAL START-UP

Before the initial start-up, remove all locking means serving for fixing the apparatus during the transport.

First actuation can be done only by professional specialist.

ANY MODIFICATION OF THIS EQUIPMENT IS FORBIDDEN!

**This equipment cannot be used nearby other instruments. If this equipment is used nearby other instruments, the equipment must be observed in order to verify normal operations in the configuration it will be used.
Instruments may be affected electro-magnetically!**

6.1. Environmental Conditions

- The apparatus may only be installed and operated in dry, well ventilated and dust-free rooms.
- The mobile dental aspirator should be installed with regard to an easy access for the operation and maintenance and good accessible rating plate.
- The apparatus should stand on a plane, sufficiently stable base (be aware of the aspirator weight, see item 3. Technical data).

Neither the lines for the connection to the electric network, nor the suction tubings may be cranked.

- To guarantee a trouble-free operation of the aspirator, the room temperature should never drop under +5°C and exceed +40°C. The ideal ambient temperatures are between +10°C up to +25°C.
- Most of the electric energy used by the suction aggregate (9) is converted to heat and given off to the ambient. During a longer lasting work, above all with the desalivating mouthpiece, the box temperature increases over 40°C causing an automatic actuation of the cooling fan. After the room is cooled down under approximately 32°C, the fan turns off again.

6.2. Electric Connection

The dental aspirator model 230V is delivered with the safety plug. It is unavoidable to observe the rules of the local electricity works. The network voltage and frequency must be in compliance with the rating plate data.

- In case that the apparatus is fix connected to the supply of electric voltage, a disconnecting device must be situated near the apparatus.
- If the apparatus is connected to the supply of electric voltage by means of a plug, the socket must be good accessible for safety reasons, so that the apparatus can safely be unplugged from the network in case of a danger.
- The network distribution must be protected by maximum 10 A.

The power cord can only be replaced by maintenance personnel!

6.3. Initial Start-up

After unpacking, put the aspirator on the floor, open the door and check whether the separator's cover (10) with sensing electrodes is attached to the separation vessel (11). In case of need, fix the cover to the vessel and then follow the instructions in the Chapter No. 9. Close the door and connect the aspirator's power cord to the socket. Change the main switch (5) over into the position "I", the indicating lamp for the network (3) goes on and the aspirator is ready for the operation.

The tips are not provided in the sterile status! It is necessary to wash and sterilize the delivered endings before use!(see chap. 9.6)

USE

7. OPERATION

THE APPARATUS MAY ONLY BE OPERATED BY THE TRAINED STAFF!

Disconnect from the network in case of a danger (unplug) ! Hot surfaces of the suction aggregate. Do not touch - Danger of burn.

7.1. Switching on the Dental Aspirator

Switch on the dental aspirator by turning the power-supply switch (5) in the position „I“. The aspirator is put into operation automatically by taking the suction tubing (6) with a mouthpiece (1) off the holder (2). It works uninterruptedly until both tubings are put on the holders again or until the separation vessel (11) is full (during an uniform operation it is filled within about 6 – 10 hours). Suction tips holder is equipped with vacuum regulation of suction unit depends on dentist needs.

The state of filling of the separation vessel is evaluated by the separation automatics which interrupts the suction by the aggregate (9) and filling of the separation vessel is indicated automatically with the indicator lamp (4) going on. At this point it is necessary to put the tubings with suction mouthpieces into the holder and switch the power-supply switch off. Afterwards open the door by pulling the grips situated on the side walls, uncouple rubber fixing lugs from the separator cover (10) and pull the separator's vessel out. Hang the separator's cover into the holder on the left side (8).

Empty the separation vessel content into the sewer, rinse the vessel with water, connect with the separator's cover in a reversal order. Set the separator into a bowl (12) in the aspirator box (the bowl is detachable). Check the connection of the cover with the vessel, seating of the separation vessel in the bowl and close the door.

7.2. Use of antifoaming tablets

Under specific exhaust conditions, increased exhaust condensate foaming may occur due to the switching off of the exhaust unit when the separation tank is not full. Antifoaming tables must be inserted in order to prevent this from occurring (see the basic accessories); add these tables to the intake screen (7). The tablets will then gradually dissolve in the input screen and will significantly decrease exhaust condensate foaming; the tablets also function as disinfection agents.

Schedule for adding tablets: Insert 1 or 2 antifoaming tablets into the intake screen (7) at the end of work and after cleaning and disinfecting the equipment.

**The tips are not delivered in a sterile state!
Replace used suction mouthpieces on tubings after each patient!**

**The used tips that are designed for multiple use are necessary to be sterilized after each patient !
The cannulas are to be exchanged after maximum of 400 cycles of sterilisation or after one year based on whichever occurs first.**

It's forbidden to cover all blow-holes placed on upper sides of unit.

Aspirator contains no spare power supply!

Before each switch on running unit it's necessary to check if door of box is closed.

8. MAINTENANCE INTERVALS - USER / TECHNICIAN

Maintenance to be done	Chapter	Time interval	To be performed by
Inlet Sieve Cleaning	9.1	per day	user
Disinfection of Tubings and Separation Vessel	9.2	per day	user
Cleaning and disinfection of the exterior surfaces of the product	9.3	at need	user
Replacement of output pre-filter	9.4	every 3 months	user
Replacement of output filter	9.5	per annum	user
Exchange of suction cannula	9.6	Maximum of 400 cycles of sterilisation or once a year *	user
Regular inspection		1x in 2 years	service technician

* Note: the suction end-pieces are exchanged at one of these events, whichever occurs first.

9. MAINTENANCE, CLEANING AND DISINFECTION

Repair works exceeding the scope of a general maintenance may only be carried out by a qualified specialist or by the manufacturer's customer service.

Use only spare parts and accessories permitted by the manufacturer. Before making any maintenance, repair or cleaning, it is unavoidable to disconnect the apparatus from the network (unplug).

Following works must be carried out from the hygiene and the correct aspirator function points of view.

9.1. Inlet Sieve Cleaning

Solid particles (sucked together with a liquid component and the air) are entrapped in the inlet sieve (7) during the apparatus operation (Enclosure No. 1, Fig. 1) which therefore must be regularly cleaned, by section 8 (though, always at the end of each shift).

The inlet sieve cleaning is possible only after the moisture inside suction tubings (6) and a sieve is reduced by means of the several seconds lasting air suction through tubings taken off the holder (2). Then turn the power-supply switch (5) off in the position „O“ and slightly lift the cover of the inlet sieve. Afterwards pull the sieve out of its place gripping the holder, remove solid particles and clean the sieve. Then put the sieve on the previous place into the holder and carry out the assembly in the contrary manner to the dismantling.

If the suction of amalgam particles is expected, it is necessary to empty the inlet sieve content into a closing vessel and forward it to the collecting centres.

9.2. Disinfection of Tubings and Separation Vessel

For the separation vessel (11) disinfection it is necessary to empty the same, rinse with warm water and clean mechanically with warm water and a disinfectant.

The disinfection of internal suction tubings and of the separation vessel should be performed by section 8, though always at the end of a shift by suction of warm water containing a disinfectant having a limited foaming power through both suction tubings (6) separately. The used disinfectant must be approved in compliance with a valid national legislation on usage for disinfected surface and a type of material. For the disinfectant use, follow the manufacturer's instructions.

9.3. Cleaning and disinfection of the exterior surfaces of the product

Use neutral agents for cleaning and disinfection of the external surfaces of the product.

The use of aggressive cleaning agents and disinfectants comprising alcohol and chlorides may lead to the damage of the surface and the discoloration of the product.

To find out whether the aspirator works correctly, the following works should be carried out in the defined maintenance intervals:

9.4. Replacement of the Output Pre-filter

It is necessary to replace the output pre-filter (15) by section 8. Before making any replacements turn the aspirator power-supply switch (5) off. The pre-filter together with the filter (14) are situated on the left lower part of the aspirator's bogie having the symbol

For its dismantling, catch the filter holder (13) with one hand and with the other hand slew the fixing lugs (16) by 90° by simultaneous pulling down. Remove the filter holder together with the output filter. Then release the output pre-filter from the cavity wherein also the output filter was situated (check also the state of pollution of the output filter).

For the reassembly, set the output filter into the filter holder, put the output pre-filter onto the filter (direct the output pre-filter with the gummed cloth towards the filter) and insert all the assembly into the cavity space from below. Slew the fixing lugs back by 90° so that they engage into the recessed parts of the filter holder.

9.5. Replacement of Output Filter

During the regular use of the device, the output filter (14) must be replaced according to the interval in Chapter 8. The procedure of disassembly and assembly of the filter is the same as for the replacement of output prefilter (15).

9.6. Sterilisation Process for Cannulas (pursuant to STN EN ISO 17664)

TYPE / MANUFACTURER OF CANNULAS DELIVERED FOR THE PRODUCT	10, 11 /CattaniS.p.A., VIA NATTA, 6/A - 43100 PHARMA – ITALY- TEL:+39 0521 607604
LIFE CYCLE OF CANNULAS DELIVERED FOR THE PRODUCT	The cannulas must be replaced after max. 400 cycles of sterilisations or after one year, whichever is sooner.
WARNINGS	Pay increased attention to the cleaning of cannulas! Do not exceed temperature of 135°C! In the case the sterilisation procedure is not observed it may come to the damage of health of a patient!
PREPARATION FOR STERILISATION	Dismount the cannulas from device – aspirator.
CLEANING: AUTOMATED	Device: Disinfection washing machine, detergent such as PULI-JET Plus (Cattani, Italy) 3. Start the cycle, wash at least 20 minutes at temperature of 60°C. 4. When taking out, check the opening of the cannula, whether the visible dirt was totally removed. If necessary, repeat the cycle or use manual cleaning.
DISINFECTION	Disinfecting agent solution such as PULI-JET Plus (Cattani, Italy) should be used in accordance with the manuals on marking. During automated cleaning the final rinse may take place at temperature of 60°C, so that thermal disinfection took place.
CLEANING: MANUAL	Equipment: Cleaning and disinfecting agent PULI-JET Plus, (Cattani, Italy), brush*, running water. Method : 4. Rinse off the surface dirt from the cannula. 5. Using a brush* apply detergent solution onto all surfaces. 6. Rinse with clean running water for 2 minutes, take care so that running water flow through the opening of the cannula. * NOTE – Using a suitable brush clean the opening of cannulas lengthwise in order to reach the entire depth of the cannula.
DRYING	If drying takes place as a part of the cycle in the disinfection washing machine, do not exceed temperature of 120°C
STERILISATION	Steam autoclave, minimally 10 minutes at temperature of 134°C. Do not exceed temperature of 135°C!
STORAGE	Cannulas should be stored in closed boxes or in suitable packages protecting them against microbial contamination after sterilisation.
ADDITIONAL INFORMATION	During the sterilisation do not exceed the maximal capacity of sterilizer and maximal temperature of sterilisation!

10. PUTTING OUT OF OPERATION

When the dental aspirator is not used for a longer time, it is necessary to carry out the cleaning and disinfection of all parts as described in the Articles No. 9.1, 9.2 and then turn the power-supply switch (5) on in the position “I”, take the suction tubings (6) off the holder (2) and let the air intake flow through (for approximately 15-20 min) so that the suction system is fully dried out. Afterwards put the tubings into the holder and turn off the power-supply switch in the position “0”, disconnect the apparatus from the network (unplug) and remove the suction mouthpieces (1) from suction tubings.

DISPOSAL

11. DISPOSAL OF APPLIANCE

- Disconnect the appliance from mains.
- To follow the rules of personal hygiene for works with contaminated material
- Clean device according, section 9
- To separate, label, packing and providing for decontamination of contaminated parts by course of national regulations
- Liquidate dental suction unit according to local effectives statutes.

Inside parts of unit can be contaminated with biological material by reason incorrect using. Before clearing and waste disposal pass on special institution for decontamination.

TROUBLESHOOTING

12. INFORMATION ON REPAIR SERVICE

Guarantee and post-guarantee repairs are ensured by manufacturer or organizations and repairers denoted by the supplier.

Warning !

Manufacturer reserves the right to make changes on the appliance that however will have no impact on substantial characteristics of the appliance.

13. SOLVING COMMON PROBLEMS

Disconnect the appliance from the network before making any interventions therein.

Works in connection with the removal of failures may only be carried out by a trained stomer service technician.

In case of repair of parts of device which might be contaminated please follow bellow mentioned instruction:

- Disconnect the appliance from mains
- To follow the rules of personal hygiene for works with contaminated material
- Clean device according, section 9
- To separate, label, packing and providing for decontamination of contaminated parts by course of national regulations
- To perform the repair of damaged parts

Failure	Problem and possible reasons	Way of removal
Aspirator does not work	Mains voltage absent	Check mains voltage in the socket, activate a protective element in the electric distribution (fuse, circuit breaker)
	Interrupted supply of el. energy Power-supply switch off	Check the power-supply switch position („I“), switch on the power-supply switch (network indicator goes on)
	Interrupted power cord	Replace the defective part
	Loose clamp of terminal board	Tighten the clamps
	Interrupted motor wiring, defective thermal protection	Replace suction aggregate
	Defective mains fuse of the apparatus	Check the state of apparatus mains fuse, replace defective fuse (fuses 2 x T 6.3 A are situated in the upper, narrowed part of the aspirator)
	Full separation vessel High current consumption	Empty the separation vessel content Defective starting capacitor, replace the capacitor
Aspirator is switching without taking suction tubings off	Suction unit is/was awash by liquid	Suction unit pass on repair to service – keep all rools and norms for work with contamination material
	Defective control automatics Defective microswitch in the suction tubings holder	Replace the defective part Replace the defective part
Aspirator is noisy	Box door improper closed	Close the box door
	Defective bearings of the suction aggregate	Replace the defective bearing for an identical type
	Incorrectly mounted filter holder	Mount the filter holder into a correct position
	Incorrectly mounted cover of the separation vessel Suction unit is/was overflowen by liquid.	Mount the separation vessel cover into a correct position Suction unit pass on repair to service – keep all rools and norms for work with contamination material
Aspirator efficiency reduced	Heavy polluted output pre-filter	Replace the pre-filter (see Maintenance intervals art.No.8)
	Heavy polluted output filter	Replace the filter (see Maintenance intervals art.No. 8)
	Leaks in the suction system	Check connections, seal untight connections
	Defective control valve in the suction tubings holder	Replace the defective part
	Heavy polluted inlet sieve	Clean the inlet sieve

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GUARANTEE

Certificate:

The Product is manufactured according to production documentation and approved technical conditions. The product is complete and meets all prescribed conditions.

Guarantee:

The Manufacturer provides a guarantee of the product, in duration of 24 months from the selling date, when the terms of guarantee document are followed.

- In the case of a claim for repairs under guarantee, a correctly completed guarantee document and certificate of quality and completeness of product must be submitted together with the product itself.
- During the guarantee period all failures that occur to the product due to production faults will be repaired free of charge.
- The guarantee does not apply to the product damaged during transportation or during incorrect storage.
- The guarantee does not apply to failures due to incorrect operation of the device.
- The guarantee does not apply to failures due to the use of the product for any other purposes other than which are described in the directions for use or as those which have been agreed with the manufacturer.
- The guarantee does not apply to the product when any unauthorised modifications have been made.

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